

# MICROCNX® ULT SERIES CONNECTORS

The **MicroCNX® ULT aseptic connectors** is specifically designed for the challenging conditions of biological media transfer in cell therapy and gene therapy applications. It can withstand ultra-low temperatures down to -190°C (cryogenic freezing applications) and fits directly into freezing cassettes used in CGT processing. Like the original MicroCNX standard connector, the ULT provides a simple, efficient method of connecting tubing for small-format biomanufacturing assemblies. The MicroCNX ULT comes in 1/8", 3/32", and 1/16" size connections to PVC tubing.



## SPECIFICATIONS

### OPERATING PRESSURE

Up to 43.5 psi, 3 bar

### OPERATING TEMPERATURE

34°F to 104°F (1°C to 40°C)

### STORAGE TEMPERATURE

-310°F to 140°F (-190°C to 60°C)

**WARNING:** Do not submerge connectors in liquid nitrogen. MicroCNX ULT has been tested and validated for cryogenic freezing use in gaseous phase of liquid nitrogen.

### TERMINATIONS

1/16", 3/32", 1/8" ID hose barb (1.6mm, 2.4mm, 3.2mm)

### MATERIALS

**Main Components:** Polycarbonate (white), Polyphenylsulfone (off white)

**Seals:** Silicone (clear), platinum-cured

**Protective Cover:** Polycarbonate (gold)

**Membrane:** Hydrophobic Polyethersulfone

### STERILIZATION

**Gamma:** Up to 50kGy irradiation.

**Autoclave:** One cycle up to 250°F (121°C) for 30 minutes

### FEATURES

PULL-CLICK-PULL →

Easy-to-Use →

Ultra-Low Temperature and Chemical compatibility →

Tubing Compatibility →

Low profile →

### BENEFITS

Easy to use three-step connection process reduces risk of operator error

No hassling with external machines or equipment, or mismatched tubing

Ability to be frozen down to vaporized liquid nitrogen temp (-190°C), and increased compatibility to harsh chemicals

Can be used with PVC, TPE and silicone tubing

Compact size to fit in freezing cassettes

## TYPICAL FLOW RATE

**Cv Value Range:** 0.04-0.27

Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

## DID YOU KNOW

Termination and flow path are made from PPSU to allow for increased chemical compatibility to handle harsh chemicals like liquid nitrogen and Dimethyl sulfoxide (DMSO).

Scan code to visit webpage



cpcworldwide.com/MicroCNX



**COLDER PRODUCTS COMPANY**  
U.S.A.  
PHONE: +1 (651) 645-0091  
TOLL FREE: +1 (800) 519-7633  
E-MAIL: info@cpcworldwide.com

**COLDER PRODUCTS COMPANY GMBH**  
Germany  
E-MAIL: cpcgmbh@cpcworldwide.com

**DOVER (SHANGHAI) INDUSTRIAL CO., LTD**  
ShangHai, China  
PHONE: +86 21 2411 2666  
TOLL FREE: +86 400 990 1978  
E-MAIL: asiapacific@cpcworldwide.com


# MICROCNX® SERIES DIMENSIONS

POLYPHENYLSULFONE with gold cover For autoclave or gamma irradiation applications.



TERMINATION	METRIC EQ.	PART NO.	A	B	C
1/16" HOSE BARB	1.6 mm	CNX17101HT	0.62" (15.7 mm)	2.35" (59.7 mm)	1.70" (43.2 mm)
3/32" HOSE BARB	2.4 mm	CNX17103HT	0.62" (15.7 mm)	2.46" (62.5 mm)	1.70" (43.2 mm)
1/8" HOSE BARB	3.2 mm	CNX17102HT	0.62" (15.7 mm)	2.56" (65.0 mm)	1.70" (43.2 mm)

## PRODUCT DIMENSIONS



## MICROCNX® ULT CONNECTORS ASSEMBLY PROCEDURE

### PINCH



Remove the protective cover on each half.

### CLICK



Join halves and click together.

### PULL



Pull membrane strips directly away from the connector.

WARRANTY: All sales are subject to Colder Products Company's limited express warranty set forth in the CPC catalog. Contact your local distributor or CPC Customer Service for warranty provisions.

Warning: Due to the wide variety of possible fluid media and operating conditions, unintended consequences may result from the use of this product, all of which are beyond the control of CPC. It is the user's responsibility to carefully determine and test for compatibility for use with their application. All such risks shall be assumed by the buyer.

COPYRIGHT © 2025 BY COLDER PRODUCTS COMPANY.

CPC, Colder Products Company, and Colder Products are registered trademarks with the United States Patent and Trademark Office.

For detailed trademark information, please visit: <https://www.cpcworldwide.com/Trademarks>